

A Mass Spectrometry Proteomics-Based Approach to Identify Target Receptors for Novel Extracellular tRNA Synthetase Fragments

Blythe C. Dillingham¹, Jennifer Brasseit², Björn Hegemann², Ann L. Menefee¹, Justin Rahman¹, Zhiwen Xu¹, Paul Helbling³, Leslie A. Nangle¹, Ryan A. Adams^{1*}

1. aTyr Pharma, 2. CSL Behring, 3. Dualsystems Biotech *Contact: radams@atyrpharma.com


Overview

- Purpose:** To identify target receptors and determine the biological functions of novel extracellular tRNA synthetase fragments with links to immune modulation
- Methods:** tRNA synthetase fragments were expressed recombinantly, and their binding to various human cell lines was assessed by flow cytometry. Receptor screening was completed in three cell-lines using the ligand-receptor capture technology LRC-TriCEPS followed by mass spectrometry proteomics analysis^{1,2}. Target validation was completed by siRNA knock-down and flow cytometry, and biological function was determined using a FRET-based enzyme-inhibitor assay.
- Results:** Utilizing this workflow, we successfully identified target cell-surface receptors for tRNA synthetase fragments and have gained insight into their previously unknown biological functions. In doing so, we have also created a novel approach which can be applied more broadly to identify receptor targets of extracellular proteins in an endogenous system.

Introduction


- While canonically known for their intracellular role in protein synthesis, full-length and splice or proteolytic variants of tRNA synthetases have been found to exist in the extracellular space where they may play an immunomodulatory role.
- Full-length Histidyl-tRNA synthetase (HARS) has been established as a molecule present in circulation that modulates T-cell activity, and a HARS variant has been shown to bind to Neuropilin-2 and to inhibit proinflammatory chemokines and cytokines.
- Alanyl-tRNA Synthetase (AARS) and Aspartyl-tRNA Synthetase (DARS) are also present extracellularly and have links to immune modulation; however, their receptor targets and downstream biological function remain unknown:
 - Auto-antibodies targeting AARS and other synthetases are present in rare anti-synthetase syndromes associated with inflammatory phenotypes such as myositis and interstitial lung disease³.
 - Full-length DARS protein and a DARS fragment are secreted from THP-1 Macrophages when stimulated with LPS (shown to the right).

DARS fragment secreted upon stimulation of THP-1 macrophages with LPS:




Materials and Methods

Recombinant Expression of AARS and DARS Fragments



- V5-tagged extracellular tRNA synthetase fragments AARS-1 and DARS-1 were expressed recombinantly in Expi293 or ExpiCHO cells and purified for downstream *in vitro* assays.
- Shown to the left are SDS PAGE (blue) and anti-V5 Western blots (grey) detecting purified AARS-1 and DARS-1 under non-reduced conditions.

FACS Cell Binding Assay




- AARS-1 and DARS-1 binding to unknown cell surface proteins on various cell types was tested by flow cytometry.
- Cells were incubated with 300 nM protein for one hour, followed by a 30-minute incubation with a fluorescently-conjugated α -V5 detection antibody.
- Binding was quantified as median fluorescent intensity (MFI), and fold changes were calculated relative to α -V5 detection antibody alone.

LRC-TriCEPS Mass Spectrometry Proteomics Receptor Screen

TriCEPS 3.0 Reagent

- Ligand of interest is coupled to TriCEPS reagent.
- Upon ligand binding, TriCEPS cross-links to the receptor targets via the hydrazide functional group.
- Ligand-Receptor complexes are isolated using the azide group for purification.

LRC-TriCEPS Receptor Screen



References


- Frei, Andreas P et al. "Direct identification of ligand-receptor interactions on living cells and tissues." *Nature Biotechnology* vol. 30,10 (2012): 997-1001. doi:10.1038/nbt.2354
- Sobotzki, Nadine et al. "HATRIC-based identification of receptors for orphan ligands" *Nat Commun* 9, 1519 (2018). https://doi.org/10.1038/s41467-018-03936-z
- Witt, Leah J et al. "The Diagnosis and Treatment of Antisynthetase Syndrome." *Clinical Pulmonary Medicine* vol. 23,5 (2016): 218-226. doi:10.1097/CPM.0000000000000171

Results

AARS-1 and DARS-1 Binding to Human Cells


Human cell type	Differentiation state	AARS-1	DARS-1
Monocyte THP-1	Naïve	+	+
Primary monocytes (classical)	Naïve	-	-
Primary monocytes (classical)	Activated (PMA)	++	++
Monocyte THP-1	M0 (PMA)	+	+
Monocyte THP-1	M1 (PMA/LPS/IFN γ)	-	-
Primary macrophages	All	-	+
Natural Killer NK-92	Naïve	-	++
Primary NK cells	Naïve	+	+
T cell Jurkat	Naïve	+	+
T cell Jurkat	Activated (α CD3/ α CD28)	+	+
T cell Jurkat	Activated (PMA)	-	-
Primary CD4+ T-cells	CD4+	-	-
Primary CD8+ T-cells	CD8+	+	+
Primary NK-T cells	Naïve	+	+
Glioblastoma U251	Naïve	++	++
Glioblastoma U87	Naïve	+	+
Lung adenocarcinoma A549	Naïve	++	++

Fold change: ≥ 20 , ≥ 3 , ≤ 3



- AARS-1 and DARS-1 binding, as well as a V5-tag control, to a panel of primary and immortalized cell-lines was assessed by flow cytometry – data overview shown in table to left.
- Three cell-lines were selected based on cell binding data for downstream receptor screening experiments – shown in graphs above (mean \pm SEM; n=2-3).


Receptor Candidates Identified from 3 Receptor Screens




- Multiple undisclosed candidate receptors were identified for both AARS-1 and DARS-1 (as indicated by numbers above).
- Some candidates appear to be cell-type specific, while others were identified across multiple cell-types.
- The candidates identified indicate potential involvement in multiple biological pathways including TGF- β and IGF signaling, oncology and fibrosis.
- TIMP-1 protein, a regulator of extracellular matrix (ECM) degradation, was identified as a candidate for DARS-1 and selected for follow-up experiments.

DARS-1 candidate TIMP-1: Binding validation and Functional Activity

DARS-1 binding to THP-1 M2 Macrophages is TIMP-1 dependent




DARS-1 enhances TIMP-1-mediated inhibition of MMP-1



- TIMP-1 regulates ECM degradation by inhibiting matrix metalloproteinases (MMPs).
- FRET-based MMP-1 enzyme inhibitor assay: TIMP-1 inhibits MMP-1 (IC_{50} = 10.2 μ M).
- DARS-1 enhances TIMP-1-mediated inhibition of MMP-1 activity in a dose-dependent manner (shown to the right).

Conclusions and Workflow

- Extracellular AARS and DARS fragments, AARS-1 and DARS-1, were successfully expressed recombinantly, and their target receptors were identified in endogenous systems using the ligand-receptor capture technology LRC-TriCEPS and mass spectrometry proteomics.
- The receptors identified provide new insight into the biological activity of extracellular tRNA synthetases.
- DARS-1 was shown to enhance TIMP-1-mediated inhibition of MMP-1 activity, indicating it may have a regulatory role in extracellular matrix remodeling.



Acknowledgements

This project was supported by aTyr Pharma and CSL Behring. In addition to the listed authors, we would like to acknowledge Esther Chong and Lauren Guy from aTyr Pharma for assistance with cloning and receptor screening experiments, respectively.

